

Pyrrolidinedithiocarbamate ammonium(PDTC)

Catalog No: #S3147


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pyrrolidinedithiocarbamate ammonium(PDTC)
Brief Description	Inhibitors
Purification	98.00%
Target Name	NF-κB inhibitor
Calculated MW	164.29
Formulation	C5H12N2S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:[NH4+].[S-]C(=S)N1CCCC1Pathways:NF-κBReceptor:NF-κBBoling pt:Melting pt:Solubility:DMSO:
42mg/mLAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a
while.Stock solution can be stored below -20 °C for several month.

References

1. Qin JD, et al. Pharm Biol. 2014 Nov;52(11):1460-1466.

Note: This product is for in vitro research use only