

## Salvianolic Acid C

Catalog No: #S3149


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Salvianolic Acid C
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	492.43
Formulation	C <sub>26</sub> H <sub>20</sub> O <sub>10</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C1=CC(=C(C=C1C[C@H](C(=O)O)OC(=O)/C=C/C2=C3C=C(OC3=C(C=C2)O)C4=CC(=C(C=C4)O)O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: Soluble Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Song J, et al. Biomed Chromatogr. 2016 Mar;30(3):376-83.

Note: This product is for in vitro research use only