

## Caudatin

Catalog No: #S3150


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Caudatin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	490.63
Formulation	C28H42O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC(C)/C(=C/C(=O)O[C@@H]1C[C@@H]2[C@]3(CC[C@@H](CC3=CC[C@]2([C@@]4([C@]1(CC4)(C(=O)C)O)C)O)O)C)/CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Fu XY, et al. Cell Mol Neurobiol. 2015 Oct;35(7):953-9.

Note: This product is for in vitro research use only