

## Brassinolide

Catalog No: #S3151


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Brassinolide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	480.68
Formulation	C28H48O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@@H]([C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2COC(=O)[C@@H]4[C@@]3(C[C@H]([C@H](C4)O)O)C)C)[C@H]([C@@H]([C@@H](C)C(C)O)O)OPathways:OthersReceptor:plant growth modulatorBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Jin H, et al. Phytochemistry. 2014 May;101:40-51.

Note: This product is for in vitro research use only