

Harmol

Catalog No: #S3152


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Harmol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	198.22
Formulation	C ₁₂ H ₁₀ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1=C2NC3=CC(=O)C=CC3=C2C=CN1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Gray Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Abe A, et al. Biol Pharm Bull. 2011;34(8):1264-72.

Note: This product is for in vitro research use only