

Guaiazulene

Catalog No: #S3159


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Guaiazulene
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	198.3
Formulation	C15H18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(C)C1=CC2=C(C)C=CC2=C(C)C=C1Pathways:OthersReceptor:OthersBoiling pt:153 CMelting pt:27-29?CSolubility:Soluble in PhenixinAppearance:BlackPowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Togar B, et al. J Intercult Ethnopharmacol. 2015 Jan-Mar;4(1):29-33.

Note: This product is for in vitro research use only