

## Metoprolol Succinate (ToprolX)

Catalog No: #S3165


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Metoprolol Succinate (ToprolX)
Brief Description	Inhibitors
Purification	99.70%
Target Name	Adrenergic Receptor antagonist
Calculated MW	652.82
Formulation	C34H56N2O10
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CC(C)NCC(COC1=CC=C(C=C1)CCOC)O.CC(C)NCC(COC1=CC=C(C=C1)CCOC)O.C(CC(=O)O)C(=O)O Pathways: GPCR/GProtein Receptor:  $\beta$ 1-adrenergic receptor Boiling pt: Melting pt: Solubility: Water: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Benfield P, et al. Drugs. 1986 May;31(5):376-429.

Note: This product is for in vitro research use only