

## Troglitazone

Catalog No: #S3170


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Troglitazone
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR agonist
Calculated MW	441.54
Formulation	C24H27NO5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: CC1=C(O)C(C)=C2CCC(C)(COC3=CC=C(CC4SC(=O)NC4=O)C=C3)OC2=C1 Pathways: Metabolism Receptor: PPAR $\gamma$  Boling  
 pt: Melting pt: Solubility: DMSO: 25mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yamazaki H, et al. In vitro inhibitory effects of troglitazone and its metabolites on drug oxidation activities of human cytochrome P450 enzymes: comparison with pioglitazone and rosiglitazone.

Note: This product is for in vitro research use only