


Description

Product Name	GSK137647A
Brief Description	Inhibitors
Purification	99.26%
Target Name	GPR agonist
Calculated MW	305.39
Formulation	C16H19NO3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: COC1=CC=C(C=C1)S(=O)(=O)NC1=C(C)C=C(C)C=C1
Pathways: Endocrinology/Hormones
Receptor: GPR120
Boiling pt: Melting pt: Solubility: DMSO: 35mg/mL
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sparks SM. *Etal. Bioorg Med Chem Lett.* 2014 Jul 15;24(14):3100-3.

Note: This product is for in vitro research use only