

Trovirdine (LY-300046)

Catalog No: #S3172


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trovirdine (LY-300046)
Brief Description	Inhibitors
Purification	99.00%
Target Name	HIV Protease inhibitor
Calculated MW	337.24
Formulation	C ₁₃ H ₁₃ BrN ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: BrC1=CN=C(NC(=S)NCCC2=CC=CC=N2)C=C1 Pathways: Microbiology&virology;

Proteases/Proteasome Receptor: HIV RT Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Uckun FM, et al. Bioorg Med Chem Lett. 1999 Sep 20;9(18):2721-6.

Note: This product is for in vitro research use only