

Taprenepag

Catalog No: #S3174


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Taprenepag
Brief Description	Inhibitors
Purification	99.00%
Target Name	Prostaglandin Receptor agonist
Calculated MW	478.52
Formulation	C ₂₄ H ₂₂ N ₄ O ₅ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC(=O)COC1=CC=CC(CN(CC2=CC=C(C=C2)N2C=CC=N2)S(=O)(=O)C2=CN=CC=C2)=C1 Pathways: GPCR/G Protein Receptor: PGE2 Boling pt: Melting pt: Solubility: DMSO: 30mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Prasanna G, et al. Exp Eye Res. 2011 Sep;93(3):256-64.

Note: This product is for in vitro research use only