

L-778123 hydrochloride

Catalog No: #S3182


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-778123 hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Transferase inhibitor
Calculated MW	442.34
Formulation	C ₂₂ H ₂₁ Cl ₂ N ₅ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: Cl.C1C1=CC=CC(=C1)N1CCN(CC2=CN=CN2CC2=CC=C(C=C2)C#N)CC1=O Pathways: Metabolism Receptor: Farnesyl transferase (FTase) Boiling pt: Melting pt: Solubility: DMSO: 25mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lobell RB, et al. Mol Y Ther. 2002 Jul;1(9):747-58.

Note: This product is for in vitro research use only