

Kevetrin hydrochloride

Catalog No: #S3184


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kevetrin hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Mdm2 activator; p53 activator
Calculated MW	179.67
Formulation	C ₅ H ₁₀ ClN ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: Cl.NC(=N)SCCCC#N Pathways: Apoptosis Receptor: Mdm2; p53 Boiling pt: Melting pt: Solubility: DMSO:

40mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only