


Description

Product Name	R112
Brief Description	Inhibitors
Purification	99.00%
Target Name	Syk inhibitor
Calculated MW	312.3
Formulation	C16H13FN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC1=CC=CC(NC2=NC=C(F)C(NC3=CC(O)=CC=C3)=N2)=C1
Pathways: Angiogenesis
Receptor: Syk
Boiling pt: Melting pt: Solubility: DMSO: 40mg/mL
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rossi AB, et al. J Allergy Clin Immunol. 2006 Sep;118(3):749-55.

Note: This product is for in vitro research use only