


Description

Product Name	PIK-inhibitors
Brief Description	Inhibitors
Purification	96.98%
Target Name	PI3K inhibitor
Calculated MW	347.37058
Formulation	C19H17N5O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:N/APathways:PI3K/Akt/mTOR signalingReceptor:PI3KBoling pt:Melting pt:Solubility:DMSO: 12 mMAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Saw S,et al. Eur J Pharmacol. 2016 Apr 5;776:176-84.

Note: This product is for in vitro research use only