

## Cyantraniliprole

Catalog No: #S3191


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cyantraniliprole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	473.71
Formulation	C19H14BrClN6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: CNC(=O)C1=C(NC(=O)C2=CC(Br)=NN2C2=C(Cl)C=CC=N2)C(C)=CC(=C1)C#N Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Bielza P, Guilln J. Pest Manag Sci. 2015 Aug;71(8):1068-74.

Note: This product is for in vitro research use only