


Description

Product Name	Ro 46-2005
Brief Description	Inhibitors
Purification	98.00%
Target Name	Endothelin Receptor antagonist
Calculated MW	473.54
Formulation	C23H27N3O6S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:COC1=CC=CC(OC2=C(NS(=O)(=O)C3=CC=C(C=C3)C(C)(C)C)N=CN=C2OCCO)=C1Pathways:GPCR/G ProteinReceptor:ET-A;

ET-BBoiling pt:Melting pt:Solubility:DMSO: 30mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Breu V, et al. FEBS Lett. 1993 Nov 15;334(2):210-4.

Note: This product is for in vitro research use only