

Daprodustat

Catalog No: #S3197


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Daprodustat
Brief Description	Inhibitors
Purification	97.00%
Target Name	HIF inhibitor
Calculated MW	393.43
Formulation	C19H27N3O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OC(=O)CNC(=O)C1C(=O)N(C2CCCCC2)C(=O)N(C2CCCCC2)C1=O
 Pathways: Angiogenesis;
 Chromatin/Epigenetic Receptor: HIF-PH
 Boiling pt: Melting pt: Solubility: DMSO: 30mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Brigandi RA, et al. Am J Kidney Dis. 2016 Jun;67(6):861-71.

Note: This product is for in vitro research use only