

## PRE-084 hydrochloride

Catalog No: #S3198


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PRE-084 hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sigma receptor agonist
Calculated MW	353.88
Formulation	C19H28ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: Cl.O=C(OCCN1CCOCC1)C1(CCCCC1)C1=CC=CC=C1  
 Pathways: GPCR/G Protein Receptor:  $\sigma 1$   
 Boiling pt: Melting pt: Solubility: DMSO: 34mg/mL  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Su TP, et al. J Pharmacol Exp Ther. 1991 Nov;259(2):543-50.

Note: This product is for in vitro research use only