


Description

Product Name	DB07268
Brief Description	Inhibitors
Purification	99.00%
Target Name	JNK inhibitor
Calculated MW	321.33
Formulation	C17H15N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: NC(=O)C1=CC=CC=C1NC1=NC(NC2=CC(O)=CC=C2)=NC=C1
 Pathways: MAPK Receptor: JNK1
 Boling pt: Melting pt:
 Solubility: DMSO: 10 mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Liu M, et al. *Bioorg Med Chem Lett*. 2007 Feb 1;17(3):668-72.

Note: This product is for in vitro research use only