

4-Keto 13-cis-Retinoic Acid Methyl Ester

Catalog No: #S3203


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Keto 13-cis-Retinoic Acid Methyl Ester
Brief Description	Inhibitors
Purification	99.08%
Target Name	Others
Calculated MW	314.42
Formulation	C ₂₀ H ₂₆ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CC1=C(C(CCC1=O)(C)C)/C=C/C(=C/C=C/C(=C\C(=O)OC)/C)/C
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt:
 Solubility: DMSO: 10 mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only