

Histone Deacetylase Inhibitor IV

Catalog No: #S3206


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone Deacetylase Inhibitor IV
Brief Description	Inhibitors
Purification	99.00%
Target Name	HDAC inhibitor
Calculated MW	325.4
Formulation	C19H23N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NC1=CC=CC=C1NC(=O)CCCCC(=O)NC1=CC=CC=C1
Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair;
NF-kB Receptor: HDAC
Boiling pt: Melting pt: Solubility: DMSO: 13 mM
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only