

Briciclib

Catalog No: #S3207


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Briciclib
Brief Description	Inhibitors
Purification	99.00%
Target Name	CDK inhibitor
Calculated MW	474.42
Formulation	C19H23O10PS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: COC1=CC=C(C=C1)CS(=O)(=O)/C=C/C2=C(C=C(C=C2OC)OC)OP(=O)(O)O
 Pathways: Cell Cycle/Checkpoint Receptor: Cyclin D1
 Boling pt: Melting pt: Solubility: DMSO: 30mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Neel Jasani, et al. DOI: 10.1158/1538-7445.

Note: This product is for in vitro research use only