

Brevianamide F


Catalog No: #S3208

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Brevianamide F
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	283.32
Formulation	C16H17N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1C[C@H]2C(=O)N[C@H](C(=O)N2C1)CC3=CNC4=CC=CC=C43
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt:
 Solubility: DMSO: 14 mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Borthwick AD. Chemical Reviews. 2012. 112 (7): 3641C3716.

Note: This product is for in vitro research use only