

PKC-412(Midostaurin)

Catalog No: #S3211


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PKC-412(Midostaurin)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor; PKC inhibitor
Calculated MW	570.64
Formulation	C35H30N4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@@]12[C@@H]([C@@H](C[C@@H](O1)N3C4=CC=CC=C4C5=C6C(=C7C8=CC=CC=C8N2C7=C53)CNC6=O)N(C)C(=O)C9=CC=CC=C9)OCPATHWAYS:Cytoskeletal Signaling; OthersReceptor:PPK; PKCα; PKCβ1; PKCβ2; PKCγBoling pt:Melting pt:Solubility:DMSO: 175.24 mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Fabbro D, et al. Pharmacol Ther, 1999, 82(2-3), 293-301.

Note: This product is for in vitro research use only