

Selexipag(NS-304)

Catalog No: #S3216


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Selexipag(NS-304)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others agonist
Calculated MW	496.62
Formulation	C26H32N4O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: OtherSMILES: CC(C)N(CCCCOCC(=O)NS(C)(=O)=O)C1=CN=C(C2=CC=CC=C2)C(=N1)C1=CC=CC=C1 Pathways: OthersReceptor: Prostacyclin receptorBoiling pt: Melting pt: Solubility: DMSO: > 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kuwano K, et al. J Pharmacol Exp Ther. 2007 Sep;322(3):1181-8.

Note: This product is for in vitro research use only