

16,17-Epoxyprogrenolone acetate

Catalog No: #S3219


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	16,17-Epoxyprogrenolone acetate
Brief Description	Inhibitors
Purification	99.13%
Target Name	Others
Calculated MW	372.5
Formulation	C23H32O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC(=O)OC[C@]12CC[C@@H](O1)C[C@@H]3[C@@]2(CC[C@H]4[C@H]3CC=C5[C@@]4(CC[C@@H](C5)OC(=O)C)C)C Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhu X, et al. Biomed Res Int. 2016;2016:5727631.

Note: This product is for in vitro research use only