

## 2-(4-(4-methoxybenzyloxy)phenyl)acetonitrile

Catalog No: #S3221


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2-(4-(4-methoxybenzyloxy)phenyl)acetonitrile
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	241.26
Formulation	C15H12FNO
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:FC1=CC=C(COC2=CC=C(CC#N)C=C2)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:PinkPowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only