

16,17-Epoxyprogrenol

Catalog No: #S3222


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	16,17-Epoxyprogrenol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	330.46
Formulation	C21H30O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC(=O)[C@]12[C@H](O1)C[C@@H]3[C@@]2(CC[C@H]4[C@H]3CC=C5[C@@]4(CC[C@@H](C5)O)C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only