

## 16beta-Methyl-16alpha,17alpha-epoxypregnenolone

Catalog No: #S3223


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	16beta-Methyl-16alpha,17alpha-epoxypregnenolone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	344.5
Formulation	C22H32O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CC(=O)C12C3(CCC4C(C3CC1(O2)C)CC=C5C4(CCC(C5)O)C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only