

16,17-Epoxy-21-acetoxypregnenolone

Catalog No: #S3225


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	16,17-Epoxy-21-acetoxypregnenolone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	388.5
Formulation	C23H32O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC(=O)OCC(=O)[C@]12C(O1)C[C@@H]3[C@@]2(CC[C@H]4[C@H]3CC=C5[C@@]4(CCC(C5)O)C)2 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Chloroform Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only