

## Paritaprevir

Catalog No: #S3226


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Paritaprevir
Brief Description	Inhibitors
Purification	99.66%
Target Name	HCV Protease inhibitor
Calculated MW	765.89
Formulation	C40H43N7O7S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism/

Immunology SMILES: Cc1ncc(nc1)C(=O)N[C@H]5CCCCC=C[C@@H]8C[C@]8(NC(=O)[C@@H]6C[C@@H](Oc3nc4cccc4c2cccc23)CN6C5=O)C(=O)NS(=O)(=O)C7CC7 Pathways: Microbiology&virology; Proteases/Proteasome Receptor: HCV NS3/4A Binding pt: Melting pt: Solubility: Soluble in DMSO Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Todt D, et al. J Antimicrob Chemother. 2017 Jan 18.

Note: This product is for in vitro research use only