


Description

Product Name	Danshensu
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	198.17
Formulation	C9H10O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC(=C(C=C1C[C@H](C(=O)O)O)O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Water: 6.2 mg/mL (warm)
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. YG Cao, et al. Br J Pharmacol. 2009 Jun; 157(3): 482C490.

Note: This product is for in vitro research use only