

NLRP3 Inflammasome Inhibitor I

Catalog No: #S3230


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NLRP3 Inflammasome Inhibitor I
Brief Description	Inhibitors
Purification	98.00%
Target Name	NOD inhibitor
Calculated MW	368.83
Formulation	C16H17ClN2O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=CC=C(C=C1Cl)C(=O)NCCC2=CC=C(C=C2)S(=O)(=O)N
Pathways: Immunology/Inflammation;
NF-κB Receptor: NLRP3
Boiling pt: 209-214 °C
Melting pt: 209-214 °C
Solubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only