

Solanesol

Catalog No: #S3236


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solanesol
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	631.07
Formulation	C45H74O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC(=CCC/C=C/CC/C=C/CC/C=C/CC/C=C/CC/C=C/CC/C=C/CC/C=C/CC/C=C/CO)/C)/C)/C)/C)/C)/C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:33CSolubility:Soluble in ChloroformAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Yao X, et al. Solanesol induces the expression of heme oxygenase-1 via p38 and Akt and suppresses the production of proinflammatory cytokines in RAW264.7 cells. Food Funct. 2017 Jan 25;8(1):132-141.

Note: This product is for in vitro research use only