

L- $\alpha$ -Phosphatidylcholine

Catalog No: #S3237


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	L- $\alpha$ -Phosphatidylcholine
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	758.06
Formulation	C42H80NO8P
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous

systemSMILES:CCCCCCCCCCCCCCCC(=O)OCC(COP([O-])(=O)OCC[N+](C)(C)C)OC(=O)CCCCCCCC=CCC=CCCCCPathways:OthersReceptor:OthersBoiling pt:Melting pt:57°CSolubility:Ethanol: >5mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mansour M, et al. Novel polyglycerol-dioleate based cubosomal dispersion with tailored physical characteristics for controlled delivery of ondansetron. Colloids Surf B Biointerfaces. 2017 Apr 26;156:44-54.

Note: This product is for in vitro research use only