

Lanolin


Catalog No: #S3238

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lanolin
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	N/A
Formulation	N/A
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCCCCCCCCC(=O)(OCCCCCCCCCCCCCCCC)Pathways:OthersReceptor:OthersBoiling pt:Melting pt:38-40 CSolubility:Chloroform: 0.1g/mL, clear to faintly turbid (

References

1. Sharma G, et al. Lanolin-based organogel of salicylic acid: evidences of better dermatokinetic profile in imiquimod-induced keratolytic therapy in BALB/c mice model. Drug Deliv Transl Res. 2017 Feb 21.

Note: This product is for in vitro research use only