

## Geranyl?Tiglate

Catalog No: #S3239


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |  |
|-------------------|--|
| Product Name      | Geranyl?Tiglate                                |
| Brief Description | Inhibitors                                     |
| Purification      | 95.00%   |
| Target Name       | Others   |
| Calculated MW     | 236.355  |
| Formulation       | C15H24O2                                       |
| Storage           | 3 years -20°C powder;2 years -80°C in solvent; |

## Images



## Product Description

Research Area:OtherSMILES:C/C=C(\C)/C(=O)OC/C=C(\C)/CCC=C(C)CPathways:OthersReceptor:OthersBoiling pt:197?CMelting pt:Solubility:soluble in alcohol , plant oilAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Sienkiewicz M, et al. The antibacterial activity of geranium oil against Gram-negative bacteria isolated from difficult-to-heal wounds. Burns. 2014 Aug;40(5):1046-51.

Note: This product is for in vitro research use only