

## Citronellol


Catalog No: #S3240

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Citronellol
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	156.27
Formulation	C10H20O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(CCO)CCC=C(C)C  
Pathways: Others  
Receptor: Others  
Boiling pt: 225°C  
Melting pt: 77-83°C  
Solubility: Soluble in alcohol  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Taylor WG, Schreck CE. (1985). J Pharm Sci. 74 (5): 534C539.

Note: This product is for in vitro research use only