

## N-Sulfo-glucosamine sodium salt

Catalog No: #S3244


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	N-Sulfo-glucosamine sodium salt
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	281.22
Formulation	C6H12NNaO8S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: C([C@H]([C@H]([C@@H]([C@H](C=O)NS(=O)(=O)[O-])O)O)O).[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Water Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only