

# Tabersonine


Catalog No: #S3246

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tabersonine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	336.43
Formulation	C <sub>21</sub> H <sub>24</sub> N <sub>2</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CC[C@]12CC(=C3[C@@]4([C@H]1N(CC4)CC=C2)C5=CC=CC=C5N3)C(=O)OC  
 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Chloroform Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only