

ALPHA-BOSWELLIC ACID

Catalog No: #S3248


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ALPHA-BOSWELLIC ACID
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	456.7
Formulation	C30H48O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@@]12CC[C@@]3(C=CC[C@H]4[C@@]3(CC[C@@H]5[C@@]4(CC[C@H]([C@]5(C)C(=O)O)O)C)C)[C@@H]1CC(CC2)(C)C)C

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Fathi E, et al. The Effects of Alpha Boswellic Acid on Reelin Expression and Tau Phosphorylation in Human Astrocytes. Neuromolecular Med. 2017 Mar;19(1):136-146.

Note: This product is for in vitro research use only