

Castor oil

Catalog No: #S3253


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Castor oil
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	933.45
Formulation	C57H104O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CCCCCCC(O)C/C=C\C\CCCCCCCC(=O)OCC(OC(=O)CCCCCCC/C=C\CC(O)CCCCC)COC(=O)CCCCCCC/C=C\CC(O)CCCCC
 CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:Colorless oilRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Thomas. Ullmann's Encyclopedia of Industrial Chemistry, 2005

Note: This product is for in vitro research use only