

## MonoMethyl auristatin F (MMAF)

Catalog No: #S3256


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	MonoMethyl auristatin F (MMAF)
Brief Description	Inhibitors
Purification	98.67%
Target Name	Others inhibitor
Calculated MW	731.98
Formulation	C39H65N5O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:CancerSMILES:CC[C@H](C)[C@@H]([C@@H](C(=O)N1CCC[C@H]1[C@@H]([C@@H](C)C(=O)N[C@@H](CC2=CC=CC=C2)C(=O)O)OC)OC)N(C)C(=O)[C@H](C(C)C)NC(=O)[C@H](C(C)C)NCPathways:OthersReceptor:Cytotoxicity( Karpas 299 cell)Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tai, Y. T.Blood. 123 (20): 3128C38.

Note: This product is for in vitro research use only