

Ponesimod

Catalog No: #S3258


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ponesimod
Brief Description	Inhibitors
Purification	98.01%
Target Name	S1P Receptor modulator
Calculated MW	460.97
Formulation	C23H25ClN2O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CCCN=C1N(C(=O)/C(=C/C2=CC(=C(C=C2)OC[C@@H](CO)O)Cl)/S1)C3=CC=CC=C3
 Pathways: GPCR/G
 Protein Receptor: S1P1
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Piali L, et al. J Pharmacol Exp Ther. 2011 May;337(2):547-56.

Note: This product is for in vitro research use only