

Creatine phosphate disodium salt

Catalog No: #S3262


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Creatine phosphate disodium salt
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	251.111
Formulation	C ₄ H ₈ N ₃ Na ₂ O ₅ P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CN(CC(=O)O)/C(=N/P(=O)([O-])[O-])/N.[Na+].[Na+]Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in WaterAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Schlattner, U et al. Biochimica et Biophysica Acta (BBA) - Molecular Basis of Disease. 1762 (2): 164C180.

Note: This product is for in vitro research use only