

## IK Polyclonal Antibody

Catalog No: #41054

Package Size: #41054-1 50ul #41054-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	IK Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	IK Polyclonal Antibody detects endogenous levels of IK protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human IK.
Target Name	IK
Other Names	IK; RED; RER; Protein Red; Cytokine IK; IK factor; Protein RER
Accession No.	Swiss-Prot: Q13123NCBI Gene ID: 3550
Uniprot	Q13123
GeneID	3550;
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

## Images



Western Blot analysis of HeLa cells using IK Polyclonal Antibody

---

Note: This product is for in vitro research use only