

Etilefrine hydrochloride

Catalog No: #S3264


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Etilefrine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	217.69
Formulation	C10H16ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: Cl.CCNCC(O)C1=CC(O)=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 122
CSolubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only