

## Potassium guaiacolsulfonate

Catalog No: #S3267


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Potassium guaiacolsulfonate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	502.6
Formulation	C14H16K2O11S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology/ Respiratory

system SMILES: O.[K+].[K+].COC1=C(O)C=CC(=C1)S([O-])(=O)=O.COC1=C(O)C=CC(=C1)S([O-])(=O)=O  
 Pathways: Microbiology & virology  
 Receptor: Antibiotic  
 Boiling pt: Melting pt: Solubility: Soluble in Water  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

N/A

Note: This product is for in vitro research use only