

Metadoxine

Catalog No: #S3269


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metadoxine
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor antagonist
Calculated MW	298.293
Formulation	C8H11NO3C5H7NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism/ Nervous systemSMILES:CC1=NC=C(C(=C1O)CO)CO.C1CC(=O)N[C@@H]1C(=O)OPathways:GPCR/G Protein;
NeuroscienceReceptor:5-HT2BBoling pt:Melting pt:102-104CSolubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Martinez, Diaz. Journal of International Medical Research. 30: 44-51.

Note: This product is for in vitro research use only